STORZ

KARL STORZ — ENDOSKOPE

203210 20

DUOMAT®

203210 20

DUOMAT®

203210 20

DUOMAT®

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Instruction Manual

Manual de instrucciones

203210 20

DUOMAT®

203210 20

DUOMAT®

203210 20

DUOMAT®

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**Important information for
users of KARL STORZ
Instruments**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.


This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the DUOMAT®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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**Controls, displays,
connectors, and their uses**

- ① Power switch
- ② Irrigation fitting
- ③ Suction fitting
- ④ Fitting for attaching suction line to the unit
- ⑤ Fitting for attaching irrigation line to the unit
- ⑥ Fitting for attaching irrigation line to the instrument
- ⑦ Fitting for attaching suction line to the instrument
- ⑧ Drawer for abridged instructions
- ⑨ Power supply cord receptacle
- ⑩ Power fuse holder
- ⑪ Potential equalization connector
- ⑫ Bottle stand

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
Read the instructions carefully
before operating the equipment!

Irrigation connector

Suction connector

Potential equalization connector


Applied part of type BF

Electronic information product
pollution control (China RoHS)

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Danger: Risk of explosion if used in
the presence of flammable anesthetics.

Caution: To reduce the risk of electrical
shock, do not remove cover. Refer
servicing to qualified service personnel.

Keep out of reach of patients.
Do not store liquids on or above the unit.

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Important information for users of KARL STORZ InstrumentsII	Photographs of the unitIII	Controls, displays, connectors, and their usesIV	Symbols employedV
General information3	Device description3	Safety instructions4	Warnings and cautions4
Intended use6	User qualification7	Safety precautions at the site of installation7	Safety precautions when operating the unit7
Installation and operating instructions8	Unpacking the equipment8	Basic equipment8	Installing and connecting up the unit8
Connecting suction tubing10	Connecting irrigation tubing11	Operating the unit12	Testing the unit for proper operation12
Maintenance13	Fuse replacement13	Cleaning, disinfection and sterilization14	Preliminary cleaning of reusable tubing15
Cleaning/preparation15	Manual cleaning of suction bottle16	Manual preparation of reusable tubing16	Machine preparation of reusable tubing17
Sterilization17	Steam sterilization18	Maintenance19	Servicing and repair19
Repair program20	Limitation of liability21	Manufacturer's warranty21	

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Device description

The DUOMAT® is an efficient suction/irrigation pump and, while particularly suitable for use in laparoscopic operations, it can also be employed in conjunction with KARL STORZ fiberscopes. Irrigation and suction are performed at preset constant values, thus requiring minimal operating effort. Suction/irrigation tubes with valves can be used to regulate irrigation and suction in laparoscopic applications.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

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Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: Risk of explosion if used in the presence of flammable anesthetics.

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Warnings and cautions

Warning: The electrical installation of the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

Warning: Keep out of reach of patients.

Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

Warning: Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Always unplug the unit before doing any maintenance work on it (e. g. cleaning).

Warning: To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

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Caution: The DUOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

Caution: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

Caution: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

Caution: Use only fuses of the correct rating.

Caution: To avoid contaminating the unit a bacterial filter must be fitted between the fluid collection bottle and the unit.

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Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Note: This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The DUOMAT® is a suction/irrigation pump for suction and irrigation in endoscopic interventions. It is applied either via special suction/irrigation instruments or in conjunction with fiberscopes. Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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Caution: The DUOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
User qualification

The DUOMAT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

The unit is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment with the DUOMAT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the instrument.

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Unpacking the equipment

Carefully unpack the DUOMAT® and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Basic equipment

- 1 DUOMAT®
- 1 Power cord
- 1 Instruction manual
- 1 Plastic wallet for instruction manual
- 1 Operator's quick reference guide (self-adhesive)

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
Installing and connecting up the unit

Note: The DUOMAT® may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

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Warning: The unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Installation and operating instructions

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The unit is equipped with a connector ⑪ for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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Set the unit on a flat surface. If a bottle holder is used, screw it onto the side of the unit.

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Remove power fuse holder ⑩ and withdraw voltage preselection module. Insert voltage preselection module into power fuse holder so that the mains voltage indicated in the window at the back of the power fuse holder agrees with the local voltage supply.

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Connect power cord. Insert power cord into power cord receptacle ⑨ as far as it will go.

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Warning: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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Note: The following installation steps depend on the accessories that you have for the unit.

Fit the bottle cap with connecting pieces onto the suction-line bottle and place the bottle into its holder.

Note: Make certain that the cap of the suction-line bottle is firmly seated since adequate irrigation will not otherwise be achieved.

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Insert the irrigation bottle into its holder.

Remove the cap of the irrigation bottle; attach the sterile bottle cap and secure it with the sealing clip.

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Connecting suction tubing

Attach the LUER-lock connection of the vacuum tube to the unit's connector ③.

Screw the bacterial filter into the vacuum line fitting ④ on the fluid collection bottle and attach the vacuum tube to the filter housing.

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Caution: To avoid contaminating the unit a bacterial filter must be fitted between the fluid collection bottle and the unit.

Join connecting piece (7) of the suction-line bottle to the instrument's suction connection by means of the second silicone tube.

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Connecting irrigation tubing

(when using the irrigation bottle supplied):

Connect the LUER-lock fitting of the short tube to the short nozzle (5) of the bottle cap. Connect the other end of the tube to the unit's fitting (2).

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Connect the long angled connecting piece of the bottle cap (6) to the suction/irrigation tube or the suction/irrigation handle by means of the second silicone tube.

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(when using sterile piercing bottles):

Introduce the puncture cannula into the irrigation solution bottle.

Connect the LUER-Lock fitting of the puncture cannula to the instrument's fitting (2).

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Connect the long, grooved section of the puncture cannula to the suction/irrigation tube or the suction/irrigation handle by means of the second silicone tube.

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Operating the unit

Switch on power switch ①.

The unit operates at preset constant pressure values for suction and irrigation.

Suction and irrigation can now be initiated using the valves on the suction/irrigation tube or on the suction/irrigation handle.

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Testing the unit for proper operation

Note: The unit must be checked to ensure safe functioning and proper condition prior to use on patients.

Remove the tube coming from fitting ② from bottle cap ⑤.

Switch on the unit. Place a finger on the end of the tube and check whether a flow of air forms.

Then seal the distal end of the suction tube with a finger and check whether a vacuum forms.

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Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Remove the power fuse holder ⑩ with a screwdriver or other tool.

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Caution: Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

	100/120 VAC	230/240 VAC
Power fuse	2 x T 0.5 A L	2 x T 0.2 A L

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Refit the power fuse holder ⑩.
Reconnect the power cord.
Test the unit for proper operation.

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Cleaning, disinfection and sterilization

Note: The instructions for "Cleaning, Sterilization and Care of KARL STORZ Instruments" must be followed. There the procedures for Cleaning, Disinfection and Sterilization are explained in detail.

Warning: Always pull out power plug before cleaning!

Caution: Avoid allowing fluids to enter the unit.

The exterior of the unit must be cleaned with a disposable cloth moistened with disinfectant (no solutions predominantly containing alcohol).

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Cleaning, disinfection and sterilization of reusable tubing

Warning: Reusable tubing must be cleaned, disinfected and sterilized using validated preparation procedures before initial use and subsequent reuse.

Warning: For all cleaning and disinfection procedures on contaminated instruments (tubing sets), the guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

Warning: Prior to the sterilization of reusable tubing, particularly in cases of manual preparation, it must be ensured that all chemical residue has been rinsed off or removed (see also the preparation steps).

Warning: Disposable tubing sets should not be reused.

Note: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical products (see Appendix, pages 29, 30).

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Preliminary cleaning of reusable tubing

Place the reusable tubing in a container filled with a cleaning solution (as per the manufacturer's instructions) immediately after use to prevent contaminants from drying on the surface of the tubing set.

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Water quality requirements

Tap water can be used to prepare the cleaning and disinfectant solutions. When rinsing for the last time for example to remove any chemical residue, microbiologically pure/sterile water must be used.

Caution: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper dilution and working times. Prolonged immersion can lead to material changes.

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Cleaning/preparation

Warning: For all cleaning and disinfection procedures on contaminated instruments, the local guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

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Manual cleaning suction bottle

Clean the suction fluid collection bottle and overflow protection unit thoroughly. Clean the channels with a brush, avoiding spatter. Then rinse thoroughly, repeating the cleaning procedure if necessary. Ensure that the floats of the overflow protection unit can move freely and do not stick due to residue.

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Manual preparation of reusable tubing

1. Disconnect the reusable tubing from the
connectors as necessary and then immerse it
in a cleaning solution (made up according to
the instructions provided by the
manufacturer).
2. Using a syringe, rinse out the tubing
underneath the (water) surface. Repeat the
procedure several times, if necessary.
Remove any contaminants using brushes,
sponges, soft cloths or a cotton applicator.
*Cleaning accessories are available from
KARL STORZ (see catalog).*
3. After cleaning, rinse the tubing and place it in
a suitable container with disinfectant solution,
making sure that all air bubbles have been
expelled. Ensure that the tubing is completely
covered in the solution for the duration of
exposure.
4. Cleaning brushes must be cleaned and
disinfected immediately.
5. Remove the tubing after the exposure time
and rinse with microbiologically pure/sterile
water. Ensure that all chemical residue in and
on the tubing is removed.
6. Dry the tubing with sterile compressed air.
7. After preparation, inspect the tubes for
cleanliness and damage. If soiling is still
visible, the preparation procedure must be
repeated.

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Maintenance Cleaning, disinfection and sterilization

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Machine preparation of reusable tubing

To prepare using a machine, connect the reusable tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed through. Visible, heavy soiling must be removed in preliminary cleaning prior to machine preparation.

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Warning: Check the entire tubing set for leaks after preparation. This makes it possible to detect a leak early, through which liquid could escape and enter the unit.

To check for leaks use leakage tester 13242 XL for example. For this purpose connect one end to the leakage tester and block off the other end, with your finger for example. If there is a leak in the tubing set, a pressure drop will be observed on the pressure gauge of the leakage tester.

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Sterilization

The suction fluid collection bottle and lid, including the overflow protection unit, can be either steam sterilized or gas sterilized. We recommend steam sterilizing the reusable tubing before use at 134°C (+3°C) using the fractionated prevacuum procedure.

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Warning: Successful sterilization is only possible on clean, disinfected surfaces.

Warning: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

Warning: Any deviations from the recommended parameters for sterilization should be validated by the user.

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Steam sterilization

KARL STORZ has validated the following sterilization method:
Fractionated prevacuum (steam sterilization)

Warning: Tubing sets may only be sterilized using the fractionated prevacuum procedure.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage. Load the sterilizer according to its instruction manual.

At the completion of the steam sterilization cycle, all instruments should remain untouched until adequately cooled.

Warning: Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

Fractionated prevacuum procedure

The fractionated prevacuum procedure consists of four phases.

- Conditioning phase
The conditioning phase removes air from the chamber by pulling a vacuum up to four times. Steam is then blown in.
- Sterilization phase
Sterilization takes place at a temperature of 134°C (+3°C) and a pressure (absolute) of 300 kpa (3 bar) over a period of 5 minutes (minimum time).
- Evacuation phase
The steam is released.
- Drying phase
Drying takes place under renewed application of a vacuum over a period of approx. 5–20 minutes.

Validated parameters for the fractionated prevacuum method:	
Temperature:	134°C (+3°C)
Pressure (p _{absolute}):	300 kpa (3 bar)
Exposure time:	5 minutes for all instruments

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Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3403 (English version).

Note: To ensure safe operation of the system, you should perform the following procedures every 12 months:

1. Check leakage current from chassis < 100 µA.
2. Check ground impedance < 0.1 ohms.
3. Check power consumption ≤ rated output.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Troubleshooting

Warning: Always unplug the unit before all maintenance work.

Symptom:
Complete failure of the unit.
No power output.

Possible causes:
No power from the power line.
Defective power fuse.

Remedy:
Check that there is electricity to the wall outlet.
Change fuses as described in the instruction manual. Make sure you fit the correct type of fuse.

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Symptom:
Poor suction.

Possible causes:
Leakage occurring in suction line.

Remedy:
Check all suction lines and replace any showing signs of leakage.
Check the seating of the sealing cap.

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Symptom:
No suction.

Possible causes:
Float of overflow protection unit blocking suction inlet.

Suggested remedy:
Check the fluid level in the vessel and empty if required. Clean the floats and check them for freedom of movement.

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Symptom:
Too little or no irrigation.

Possible causes:
Loose connection or leak in a tube. A tube is blocked or pinched off. Stop valve is defective.

Suggested remedy:
Check tubes and replace if necessary.
Send unit in for servicing.

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Standard compliance (for 20 3210 20)
According to IEC 60601-1:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type BF 

IEC 60601-1-2 :2001:

Please read the Electromagnetic Compatibility Information in the appendix (p. 31-42).

Directive compliance (for 20 3210 20)
According to Medical Device Directive (MDD):

This medical product belongs to Class II b

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

Note: The code number after the CE mark indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized spare parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their levels of technical training, to open or repair the instrument. Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
**KARL STORZ Irrigation, Suction
and Coagulation System**

KARL STORZ presents a new state-of-the-art system for laparoscopic irrigation, suction and coagulation. This new system offers an autoclavable handle which accepts suction/irrigation tubing and a wide variety of electrodes. The entire system is compact, ergonomically designed and easy to use.

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Typ

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
(A) = Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets

Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Plus	1, 2, 3
	Sekusept Pulver Classic	1, 2, 3
Estee Pharma GmbH	Ultrademit AF	1, 2, 3
	Ultrasseptin Aktiv	1, 2, 3
	Ultrasseptin Classic	1, 2, 3
	Afid	1, 2
Fresenius AG	Afid plus	2
	Sporcid FF	1, 2
Holifa Polska Sp. zo.o	Polsept Holifa	1, 2
Indeba Industria E	Letahdeido	2
Jose Collado S.A.	Darodur 4000 Liquido	1
	Darodur 9000	2
	Darodur Sinaldehyd 2000	2
	Aldasan 2000	1, 2
Lysoform	Almyrol	1, 2
Dr. Rosemann GmbH	Desoform	1, 2, 3
	Lysoformin 3000	1, 2
	MedDis	2
Medichem International	MediZyme	1
	Edisonite Super	1, 3
Merz Hygiene GmbH	Mucadont-IS	2, 3
	Mucadont-Zymaktiv	1, 3
	Mucocit-T neu	1, 2, 3
METASYS Medizintechnik GmbH (formerly LONZA)	ID 50	2
Orochemie	A 20	1, 3, A
Promagent AB	Wavacide	2
	Descoton Extra	1, 2
Schumacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2, 3
	Perfektan Neu	1, 2
	Perfektan TB	1, 2, 3, A
	Gigasept	2, 3
	Gigasept AF forte	1, 2
	Gigasept FF	1, 2
Schülke & Mayr GmbH	Gigasept Instru AF	1, 2, 3, A
	Gigasept Med forte	1, 2
	Gigazyme	1, 3
	Lysetol V	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1, 3
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
	neodisher SeptoClean	1, A
	neodisher SeptoPreClean	1, 2
Whiteley Industries PTY.Ltd	Aidal Plus	2
	Matrix	1
Weigert, Dr. GmbH & Co.		

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination shaft, **not** the housing/handle part.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information. Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
Weigert, Dr. GmbH & Co.	neodisher DuoClean	1, A
	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The DUOMAT® model 203210 20 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The DUOMAT® model 203210 20 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The DUOMAT® model 203210 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

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Table 201		
Guidance and manufacturer's declaration – electromagnetic emissions		
DUOMAT® model 203210 20 is intended for use in the electromagnetic environment specified below. The customer or user of the DUOMAT® model 203210 20 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	DUOMAT® model 203210 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The DUOMAT® model 203210 20 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202 Guidance and manufacturer's declaration – electromagnetic immunity			
DUOMAT® model 203210 20 is intended for use in the electromagnetic environment specified below. The customer or user of the DUOMAT® model 203210 20 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

DUOMAT® model **203210 20** is intended for use in the electromagnetic environment specified below. The customer or user of the DUOMAT® model **203210 20** should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V_{rms} 150 kHz to < 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V_{rms}</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the DUOMAT® model 203210 20, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = [3,5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3,5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than **3 V/m**.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
DUOMAT® model 203210 20**

The DUOMAT® model 203210 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:
Company stamp/signature:

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Tipo del aparato: Número de serie:

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Firma/Fecha:

To be filled out by instrument owner:
Return address/company stamp:

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KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenixia Street, Mina Elhoshn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax: +961 1 365151

KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O. Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kazakhstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51, Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 — 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com